

## PROTOCOL FOR COUPLING ANTIBODIES TO POLYSTYRENE MICROSPHERES

### Prepare in advance:

- MES buffer: 100 mM MES (2-(N-Morpholino)ethanesulfonic acid), pH=6
- Glycine buffer: 100 mM Glycine, adjust pH=6 with diluted NaOH and HCl
- Pluronic F127: 10 % Pluronic F-127 in milliQ H<sub>2</sub>O, filter in steriflip
  
- Polystyrene microspheres (1  $\mu$ m diameter), amino-modified, fluorescent:
  - SigmaAldrich L2778 (1 $\mu$ m, NH<sub>2</sub>, red fluorescent)
  - SigmaAldrich L1030 (1 $\mu$ m, NH<sub>2</sub>, yellow-green fluorescent)
  
- EDC: *N*-Ethyl-*N'*-(3-dimethylaminopropyl)carbodiimide hydrochloride (Sigma). Store tightly wrapped at -20C

1. Dilute 25  $\mu$ l microspheres in 1 ml MES buffer
2. Sonicate 2 min and centrifuge 5 min at 10,000 g
3. Resuspend in 800  $\mu$ l MES buffer and vortex well
4. Add 10-50  $\mu$ g antibody in no more than 100  $\mu$ l water or a suitable buffer without free amino groups (make sure the Ab is not in Tris buffer!!). As negative control, use 50  $\mu$ g bovine serum albumin prepared at 0.5 mg/ml in H<sub>2</sub>O.
5. Prepare fresh EDC (Dimethylaminopropyl-3-ethylcarbodiimide) only at this point!  
10 mg EDC in 1 ml MES (vortex well)
6. Add 100  $\mu$ l EDC to the beads-Ab mixture and agitate well.
7. Keep at room temperature x 2 hours with constant rotation. Agitate and sonicate briefly every 20 min.
8. Quench the reaction with 350  $\mu$ l glycine buffer, 30 min x room temperature
9. Centrifuge 5 min at 10,000 g
10. Resuspend in 1 ml PBS with 0.2 % Pluronic F127. Mix well until the suspension is homogeneous.
11. Centrifuge 5 min at 10,000 g. Repeat steps 10-11 again.
12. Resuspend in 1 ml PBS containing 0.2 % Pluronic F127 and 0.1 % BSA
13. Sonicate x 2 min.
14. Store in aliquots at 4C, wrapped in aluminum foil.
15. Use within two weeks (dilution 1/20 – 1/200 to detect surface antigens in cultured cells).

Bibliography: Riccio et. al. (1997) Science 277, 1097  
Heerssen et. al. (2004) Nature Neuroscience 7, 596